

Ward's Custom Microscope Slides

Choose Your Customization Options

Our complete range of custom slides capabilities includes a variety of tissue sampling, sectioning, mounting and staining options, as well as flexible turnaround times to meet your needs.

Tissue Sample Provided	Processing Options	Sectioning Options	Mounting Options on Glass Slide	Staining & Covering Options for Glass Mount	Requested Turnaround Time
Fixed	Embedded in wax (standard is one specimen per block)	5–20 μm standard	Single specimen per slide; Multiple specimen per slide; Serial section	Dehydrate & stain (H&E standard; other options)	3–4 weeks
Frozen in embedding medium (orientation by customer)	Embedded in plastic resin (JB101)	1–2 μm standard	Single specimen per slide	Stain (H&E standard)	4–6 weeks
	Free-floating section	50 μm (can be sent back in sucrose)	Float mount		6–8 weeks
	Cryostat	10–20 μm	Dry mount		Not time sensitive

Don't see what you need here?

Contact us for additional capabilities upon request.

Email sciencehelp@vwr.com or call 800-962-2660.